

May.2013

## **VI2010 CDS**

### **User Guide**

Tel/Fax: +86-28021919

Tel/Fax: +86-28021920

Email: sales@hplc.com.cn

[www.surwit.com/en](http://www.surwit.com/en)

VI2010 CDS User Guide

ENG

**© 2013 Surwit Technology Inc.**

## Table of Contents

1.0 Introduction.....	3
1.1 Specifications .....	3
1.1.1 System Requirements.....	3
1.1.2 CDS Main Specification .....	3
1.2 Installation.....	4
1.2.1 Installation.....	4
1.2.2 Uninstall .....	4
1.3 VI2010 Files and Directories .....	4
1.3.1 Single Needle (*.SS) .....	4
1.3.2 Instrument Configuration(*.INC).....	5
1.3.3 Spectra Processing (*.CP).....	5
1.3.4 Spectra (*.CDF file) .....	5
1.3.5 Calibration Curve(*.CAL) .....	5
1.3.6 LOG file (*.LOG) .....	5
1.3.7 Report style files (*.STY files).....	5
2.0 Spectra operate Mode.....	6
2.1 Zoom in .....	6
2.2 Zoom out .....	6
2.3 Drag: .....	6
3.0 Introduction.....	6
3.1 Five system icons in the left of Main window, User Account, Configuration, Log, .....	6
3.2 User Account .....	7
3.3 System Configuration .....	8
4.0 Audit.....	10
5.0 Report style .....	11
6.0 Data Acquisition.....	12
6.1 Data Acquisition.....	12
6.2 Spectra processing.....	15
6.3 Introduction of Menu Bar.....	16
6.4 Spectrum .....	17
6.5 Calibration.....	27
6.6 Menu Introduction.....	28
7.0 Calibration Curve Made Procession.....	28
7.1 ESTD.....	28
7.2 ISTD.....	30
8.0 Report.....	31

## **1.0 Introduction**

Thank you for buying Surwit VI2010 Chromatography Data System, in order to ensure your proper operation, please read this manual carefully, it has introduced the main performance and technical indicators, installation, unloading, debugging, and operation example, etc. If has any other questions, please contact with us for the timely answer.

Any product has shortcomings, our product is not exceptional. Please contact with us if you have any suggestion. Please contact us for improving the service.

## **1.1 Specifications**

### **1.1.1 System Requirements**

Hardware Computer

CPU: 166 MHz or faster

Memory: 32M Bytes or more

Disk Drive: 1 hard disk, 1 CD-ROM drive

USB: 2 port

Monitor: VGA display and graphics card

256 colors, 1024 by 768 pixels

Software Windows: Microsoft Windows 98/2000/XP/VISTA/WIN7 or later

Internet Explorer 4.0 or later

### **1.1.2 CDS Main Specification**

General

Number of Channels: 2

Input Voltage: -2.5V ~ 2.5 V

Input Resistance: > 10 M

Sensitivity of Integration: 1 u V/s

Minimum resolution: 1 u V

Dynamic Range:  $10^7$

Linearity:  $\pm 0.1\%$

Peak Processing

Number of peaks: > 1000

Width of peak: 0.01 s

Automatic Time Programming

Manual Integration enabled

Science Technology (Hangzhou) Inc.  
Tel: 86-0571-28021919

[www.surwit.com](http://www.surwit.com)      Email: [sales@hplc.com.cn](mailto:sales@hplc.com.cn)  
Add: 3 B. No.611 Dongguan Road, Binjiang District, Hangzhou

Automatic identification of complex peaks and precise partition of overlapped peak  
Automatic tracing and correcting base line  
Automatically eliminate the affect of negative peak  
Meet the GMP regulation

#### Methods for Identifying Peak

Conservative time method  
Components table method

#### Parameters for Integration

Peak area  
Peak height

#### Method for Quantitative Calculation

1. Normalized Method,
2. Normalized Method with Proportional Factor,
3. Internal Standard Method,
4. Grouping Method,
5. Multiple Internal Standard Method,
6. External Standard Method,
7. Logarithmic Method.

## 1.2 Installation

### 1.2.1 Installation

Before installation, please check if there is NET 2.0 Platform, if not, click the System

File in CD, double click  dotnetfx to set up.

Double click set up, choose the save path, then click “Next” to finish, if install in C, you should make folder in C.

### 1.2.2 Uninstall

Click the in folder to uninstall, please save the data to avoid lost before uninstall.

## 1.3 VI2010 Files and Directories

### 1.3.1 Single Needle (\*.SS)

Single needle is a workstation for chromatograph signal acquisition sequence. Involves saving path, sample number, condition, file name, acquisition time, equipment configuration files, spectra processing files and so on. Single needle is saved as a \* SS File, can be directly opened to use before sampling.

### **1.3.2 Instrument Configuration(\*.INC)**

Instrument configuration is the record of the experiment information, such as chromatographic column, mobile phase, flow and external control use, acquisition frequency. Counteraccusation version can also set the pump flow, gradient, pressure and the wavelength of the detector, etc. Save as \*. Inc can be used under the same equipment configuration condition.

### **1.3.3 Spectra Processing (\*.CP)**

VI2010 can export the integrals, result and performance of spectrum figure through spectrum diagram window, method can be used the CP files. In the single needle file loading spectra can be used quickly.

### **1.3.4 Spectra (\*.CDF file)**

VI2010 use the most universal CDS \*.CDF file format, AIA of other CDS file format can be also imported.

### **1.3.5 Calibration Curve(\*.CAL)**

Calibration curve which is made under the calibration window, include external standard method or internal standard method. Concentration, system suitability can be calculated by the calibration file.

### **1.3.6 LOG file (\*.LOG)**

LOG file is used to record the CDS running, such as open, close and sampling, etc. LOG file can be formed every day by the name of date.

### **1.3.7 Report style files (\*. STY files)**

VI2010 Report style can be set in the main window, all the report style is saved as \*. STY file

## **2.0 Spectra operate Mode**

### **2.1 Zoom in**

Hold the left mouse button in the spectra window from left to right, the selected area will be enlarged after let it go.

### **2.2 Zoom out**

Hold the left mouse button in the spectra window from right to left, the selected area will be the whole spectra after let it go.


### **2.3 Drag:**

Hold the right mouse button in the spectra window, the other part of spectra can be shown in the window.

## **3.0 Introduction**






### **Main Window**

Clarity Main Window include acquire, calibration, spectra, system and else buttons, click the relevant module to operate.





Input the Name and Passkey to login in.

### **3.1 Five system icons in the left of Main window, User Account, Configuration, Log, Report Style and About respectively.**

	User Account	Set visit permissions and password especially administrative user account
	Configuration	Set the configuration of the instrument.
	Audit	Shows the selected records of the mission and set their logs
	Report Style	Edit the report, set the report style.
	About	Information such as the version on CDS, ownership and else.

### 3.2 User Account

Click  , the window below will show you the dialog box.



The screenshot shows the 'User Account' dialog box. It contains several sections: 'Company' and 'Owner' fields at the top; 'User Info' section with 'Name' (Admin), 'Descript' (This is Administrator), 'Create Date' (2013年3月26日), 'Passkey Sta' [None], 'Last Modify' (0001年1月1日), 'Last Login' (2013年3月26日), and a 'Sign' area; 'Access Rights' section with checkboxes for 'User Account', 'Sys. Config.', 'Audit', 'Calib', and 'Chrom.'; and an 'Access' table with columns 'Access' and 'Ins. Name' containing two rows: '1' with 'LC' and '2' with 'LC2'. There are 'OK' and 'Cancel' buttons at the bottom.

Click the sign area, you can choose the electronic signature, support \*.JPG Format, which will be shown on the right below corner of report.


Access area can edit the different permission of different account, only authorized account can check and edit the contents.

	Access	Ins. Name
1	<input type="checkbox"/>	LC
2	<input type="checkbox"/>	LC2

### 3.3 System Configuration




Click in main window, the dialogue box will be open as below.





Connect with devices:

Counterchange will be realized after successfully connected the CDS with devices.




Double click , below window will open, click Connect, when “2239” appear, open the channel, channel lamp will be green flashing.







Before connect





After connect

Choose the needed device, click right mouse button—Equip to instru.




Ins. Num: It display the total number of instruments, click , the number will

add 1. The instruments in main window will be added as below, click  , the number will be decreased.



Instrument: Choose instrument as below, such as below :





Instrument Style can be choosed in the menu, name can be changed.  
Click instrument image, the image can be exchanged.


## **4.0 Audit**

Audit window is used for record the operation of CDS, account who has the right can

view , click  to log in.




Click , below window will be shown,



History: When the icon  appear, the file has been saved as pdf. Double click




, file can be viewed. The files are in the Document of Audit folder.

How to convert to pdf? Click , the print window will be shown, choose Microsoft XPS Document Writer, click Print, input the file name to save.

## **5.0 Report style**


Click , set up the report style in the below window,



## **6.0 Data Acquisition**

### **6.1 Data Acquisition**

Data acquisition is used to collect the signals of the instrument, which can be get in the login window.




Pump specification can be modified in the below window:

Four Unit can be choosed, include MPa , kfg/cm , bar , psi.

1MPa=10<sup>6</sup>Pa=10bar=145.03psi

1 bar=10<sup>5</sup>Pa =14.503psi

1psi=6.895kPa=0.0689476bar



Double Click the below area can changed the solvent color, the corresponding color will be shown according to the proportion.


Channel Lamp:


Blue: Open the CDS

Green: Connect to instrument







Note: if the color turned to blue in the process of running, please check the connection of CDS and PC, or check the system specification if OK from the main window.

File: click  **Options** to open the below window for setting.






Analysis:

	 Inject F3	Start to run
	 Stop F4	Stop to run
	 Abort Ctrl+Q	Quit without save

## 6.2 Spectra processing

After data acquisition, open the spectrogram, CDS will automatically identified the chromatogram peak. Click on the main window "Chrom proc" icon to enter like below:







Note: For large data spectrogram, open will cost longer time.


### 6.3 Introduction of Menu Bar

#### File

File Browser Click , file name will appears at the left side of window; double click can open the spectrogram.




Overlay Mode Choose , multiform spectrograms can be opened.

Open  \*.cdf or \*.dat ( N2000 ) style spectrogram file can be opened.

If open \*.dat ( N2000 ) file, below window will be appear: Frequency can be set up by clients.





Close 

Save 

Export \*.CDF, \*.TXT, \*.CP, \*.WMF can be exported, the files will be in the file of “Export”, the spectrogram can be output in the Word and Excel.

Apply. AP File Spectrogram can be displayed in the CP file.

Exit Exit the spectrogram window



Preview





Print



Options

## 6.4 Spectrum



Spectrum graph window size is adjustable, move the mouse to the interface of the spectra and the results table, the mouse into a , at this time hold the mouse Can adjust the size.

Hold the left mouse button drag from left to right, spectra can be amplified, from right to left shows full spectrum diagram, Hold down the right mouse button, spectra can be moved.

Integ. Integral Table display the min. peak width and the peak height.

Integ.	Result	Performance	Method	Audit	Text
	Operate	Time A	Time B	Value	
1	Global Min. Peak Width			0.010	min
2	Global Min. Threshold			0.100	mV
3	Add Positive	10.74	13.48		
4	Add Positive	1.99	2.35		
5	Add Positive	7.06	8.43		

Red color shows integral failure, does not meet the requirements for the integral. Right Click or by "Delete" on keyboard can delete the integral.

Four rows on the right can do the integral by hand.



#### Detector delay

Manually select or input lag time value, spectra will be left or right offset by the input value time range, positive to left, negative to the right.



#### Global minimum peak width

Set the minimum peak width, smaller one will not be identified.





#### Global minimum peak height

Set the minimum peak width, smaller one will not be identified.



#### Negative peak detection



#### Shoulder cut area ratio

The default value of CDS is 3, two incompletely separate peak, if their area ratio is more than 3, it was identified as shoulder peak, if less than 3 will be identified as overlapping peak.




#### Shoulder cut slope ratio

The default value is 2.5, two incompletely separate peak, if their slope is bigger than 2.5, it was identified as overlapping peak, if less than 3 will be identified as shoulder peak.




#### Add positive peak

Add positive peak within the selected scope




Add negative peak

Add negative peak within the selected scope





Seperation






According to the requirements two kinds can be choosed, automatic adjustment and unchanged baseline.





Baseline will be adjusted automatically after delete one peak.



 Baseline will not be changed




 Delete peak




 Set smallest peak width


 Set smallest peak height


 Set smallest half peak width



 Set smallest peak area


 Valley point     Add or change the original valley position within a certain range



 Make baseline across the valley point, adjust peak (group) to be single peak

 Integration baseline, adjust the peak (group) as the overlapping peak



 Adjust peak (group) the baseline level to be forward





 Adjust peak (group) the baseline level to be backward



 Before cut sets



 Tail cut



Solvent peak



Group peak

## Result

Integ.	Result	Performance	Method	Audit	Text		
%*#s							
	Compound Name	Peak Style	Retain [min]	Area [mV.s]	AreaPer [%]	Height [mV]	Amount [-]
1		Overlap	2.197	45.926	0.806	6.152	45.926
2		Overlap	2.342	14.877	0.261	1.460	14.877
3		Overlap	2.727	4.629	0.081	0.708	4.629
4		Overlap	2.865	5.720	0.100	0.407	5.720
5		Overlap	3.385	3.168	0.056	0.355	3.168
6		Overlap	3.523	1.844	0.032	0.335	1.844
7		Overlap	5.533	14.192	0.249	1.453	14.192
8		Overlap	6.148	3.837	0.067	0.131	3.837
9		Single	7.412	5457.4...	95.742	411.034	5457.4...
10		Single	9.230	2.428	0.043	0.153	2.428
11		Single	11.805	139.926	2.455	1.500	139.926
12		Single	16.943	6.150	0.108	0.107	6.150
	Total			5700.1...	100.000	423.795	5700.1...

Calib. File

Calculation ESTD

Result Style

☐ Hide Internal

☒ All Dtec. Peak

☐ All Ident. Peak

☐ All Calib. Peak

☐ Use Scale Factor

Factor

Unit

Unident Peak

Response Area


Factor  [Amt/Res]

Total Amount: ?


Volume: ?

(Comb.) Column Set			
	Vis	Column Name	Format
1	<input checked="" type="checkbox"/>	Compound Name	0.000
2	<input checked="" type="checkbox"/>	Peak Style	0.000
3	<input checked="" type="checkbox"/>	Retain[min]	0.000
4	<input type="checkbox"/>	Start[min]	0.000
5	<input type="checkbox"/>	End[min]	0.000
6	<input type="checkbox"/>	Width[min]	0.000
7	<input type="checkbox"/>	Half W.[min]	0.000
8	<input checked="" type="checkbox"/>	Area[mV.s]	0.000
9	<input checked="" type="checkbox"/>	AreaPer[%]	0.000
10	<input checked="" type="checkbox"/>	Height[mV]	0.000
11	<input type="checkbox"/>	HeightPer[%]	0.000
12	<input type="checkbox"/>	√(Area)[-]	0.000
13	<input type="checkbox"/>	√(Area)Per[%]	0.000
14	<input type="checkbox"/>	S/N[-]	0.000
15	<input checked="" type="checkbox"/>	Amount[-]	0.000

Calib. File: Calib File can be chosen in the space, which can be made in Calibration.  
Performance: All data here is calculated by the EP, 50% of peak width.



Right click can choose the items below:



**Method:**

The system specification can be choosed before collect data, please make sure all information has been confirmed.



Integ.	Result	Performance	Method	Audit	Text
<div>Info. Parameters</div> <div> <div> Descript Column Mobile Phase Flow Pressure Detect Temp. Note </div> <div> Info. Analyst [N2k.DAT] Inj. Time ? SampleID ?  Range[mV] 1.977, 413.599 Freq[Hz] 10 Method File ? </div> </div>					

**Audit:** All the operation on spectra will be recorded here to ensure the authenticity.

Integ. Result Performance Method Audit Text

User Right wouldn't be changed, after frozen

Freeze User Right

	User Name	Open	Save
1	[AnyUser]	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
2	Admin	<input type="checkbox"/>	<input type="checkbox"/>
3	Anony	<input type="checkbox"/>	<input type="checkbox"/>

X	User Name	Time	Style	Action	Note
1	Admin	2013年4月2日 15:20:16	Signal	Import	C:\Documents and Settings\Administrator\桌面\11-14 A3峰样 0009.dat

User permission can be set on the left side.

**Text:** Include all the spectra information, such as run time, end time, peak height and area and other information.

Chrom. Graph	(Comb.)Result	Summary	SST				
	Compound Name	Peak Style	Retain [min]	Area [mV.s]	AreaPer [%]	Height [mV]	Amount [-]
1		Overlap	2.197	45.926	0.806	6.152	45.926
2		Overlap	2.342	14.877	0.261	1.460	14.877
3		Overlap	2.727	4.629	0.081	0.708	4.629
4		Overlap	2.865	5.720	0.100	0.407	5.720
5		Overlap	3.385	3.168	0.056	0.355	3.168
6		Overlap	3.523	1.844	0.032	0.335	1.844
7		Overlap	5.533	14.192	0.249	1.453	14.192
8		Overlap	6.148	3.837	0.067	0.131	3.837
9		Single	7.412	5457.4...	95.742	411.034	5457.456
10		Single	9.230	2.428	0.043	0.153	2.428
11		Single	11.805	139.926	2.455	1.500	139.926
12		Single	16.943	6.150	0.108	0.107	6.150
Total				5700.1...	100.000	423.795	5700.152

**(Comb.) Result:** This table includes all the results information.

Chrom. Graph	(Comb.)Result	Summary	SST
	Chrom. Name	Amount	
1	11-14 A3标样0009		
	Mean		
	Std. Dev.		
	% RSD		

Summary: When open two or more than two spectra, input the same Chrom.Name in every spectra, or load the calibration file, the same Chrom. Average (Mean), RSD (Std. Dev.), %RSD will be displayed in the summary.

Chrom. Graph

(Comb.)Result

Summary

SST

[...]

Use	R	Compound Name	Retain [min]
-----	---	---------------	--------------

X	Chrom. Name	Retain [min]	Area [mV.s]	Height [mV]	Amount [-]
	Upper Limit :				
	Lower Limit :				
	RSD% :				
	Average				
	RSD[%]				
	Para. Result				
1	11-14 A3标样0009				

Property

☒ EP(Europ)

Descript

☐ USP(Amer)

☐ JP(Japan)

SST: To verify the chromatography system based on the chromatographic evaluation. Load the files at the open spectra, components will be refreshed in the calibration file.


### SST Making Proccession

Open the spectra, load the calibration files. Click “SST”, below window will be displayed.

Click the Use of every compound in the left, input the limit value and the RSD%, CDS will calculated and judged.

Red color means overpass the limitation, Blue color means no judgement without limitation value, others is black.

Save as \*.SST file after edited.



Without limitation value




All parameters are qualified



Parameters are more than the given limitation.








## 6.5 Calibration

Open the calibration window as below:




## 6.6 Menu Introduction


### Files

New		New calibration curve
Open		Open calibration curve
Save		Save calibration curve
Save as		Save to another files
Open Std.		Open the standard spectra which used for made calibration curve
Close Std.		Close the opened standard spectra
Preview		Preview the calibration curve
Print		Print the calibration curve
Close		Close the calibration winow

## 7.0 Calibration Curve Made Procession


### 7.1 ESTD

Open the first spectra, click add all the peaks , input the compound name,



change the Resp. Style to be Area or Height, and input the Amount correspondently.

Open another standard spectra, given the Cali. Dot. To 2 . The


first spectra can be closed or reserved. Click , input the Amount correspondently.

Note: if the Resp. Style is 0, please check if the reserve time of spectra 2 is within the spectra 1.

Then open spectra 3, given the Cali Dot To 3. , click , input the Amount correspondently.

Then open spectra 4, given the Cali Dot To 4. , click , input the Amount correspondently

Click “Compounds” in the left botton to see the calibration curve.



Menu above:

Response: Displayed the response value

Amount: Amount of compounds, which can be input in the table, the unit can be

changed by click

Calcu Value: The amount value in the calculation mode

Depart: (calculation value-amount)/amount\*100%

Recs: Number of Value which added in the same calibration point.

Use: means use this point, will be displayed in the curve, if no choice, will be displayed.

Resp. Style: Style can be choosed in this menu

Curve Style: Style can be choosed in this menu

Origin Style: Include three mode:

Ignore: ignore the point (0,0)

Join: join the point (0,0) as a point, which will be used when only one data point

Pass: pass the original point (0,0)


Remnant.: Resp.  $\frac{(\text{Resp. Value} - \text{Actual Value})}{\text{Actual Value}}$

## 7.2 ISTD

The basic step of ISTD is same with ESTD.

Notice: ISTD without equation.



Press  to save the calibration file (\*.cal).

## 8.0 Report

